



## »»Blue M

Hi-Treet Box Furnace

### Hi-Treet Box Furnace

Our Hi-Treet Box Furnace has long life and reliability with precision high-temperature uniformity. There's a minimum heating of outer skin to maintain safety while in use and also product load protection.



#### ► Features

- Standard nitrogen atmosphere, but can operate with air or argon
- Lightweight Moldatherm® insulation
- Silicon carbide hearth plate
- Side-mounted double-hinged plug door
- SSR power module
- 8" wide furnace is stackable
- Temperature
  - Up to 2500°F (1371°C)
- Chamber Sizes
  - 0.67 Cubic Feet - 4.5 Cubic Feet

Model	Chamber Dimensions (in.)			Atmos. CFH	Input kW	Voltage	Load Wt. (lbs.)	Ship Wt. (lbs.)
	Width	Depth	Height					
11-SCMT-8188-25A	8 (20.32)	18 (45.72)	8 (20.32)	50	12.5	240/1PH/60HZ	50	850
11-SCMT-182418-25A	18 (45.72)	24 (60.96)	18 (45.72)	100	40	240/3PH/60HZ	240	2,120

All specifications are subject to change without notice.

► Options:

- Furnace Stand
- Quench Tank
- Programmable Controller with Dual Digital Read-outs
- Transformer for 480V Operation



**Gruenberg, Blue M, Tenney**

Mailing Address: P.O. Box 150 | White Deer, PA 17887-0150 | USA

Phone: (570) 538-7200 | Fax: (570) 538-7380

Physical Address: 2821 Old Route 15 | New Columbia, PA 17856-9396

info@thermalproductsolutions.com

*Specifications and Product Information are subject to change without notice.*

THERMAL PRODUCT SOLUTIONS

»»Blue M »»Gruenberg »»Lindberg/MPH »»Tenney